

Technical Data Sheet

Production Name: MSEP-01

Brief Introduction

MSEP-01 is an animal-derived component-free (ADCF) cell supplement that designed for reducing the usage of fetal bovine serum (FBS). By adding to culture medium, it can substitute up to 80% FBS usage without cell adaption. It's suitable to grow different types of mammalian cell under low FBS environment over multiple passages, including myoblasts, fibroblasts, and epithelial cell. MSEP-01, extracted from edible plants, is ethical-friendly and has lower risk of infectious pathogen contamination. Currently, serum-free formulations typically have low performance across cell types and its expense is too high for scale in cultivated meat. Furthermore, MSEP-01 has synergetic effect with ITS supplement to enhance the cell proliferation, though ITS supplement is widely used to reduce the FBS consumption. MSEP-01 is a cost-saving, high performance, versatility across cell types and adaption-free cell supplement.

Application

Cell Culture Medium Supplement
Fetal Bovine Serum Reduction
Complete Medium Formulation

Instruction to Use

Reconstituting MSEP-01 with basal media intended to use. Incubate at 37 °C for 30 minutes to dissolve MSEP-01 completely. Gently pipetting to homogenize MSEP-01 solution before use. It's normal if precipitation appeared in MSEP-01 solution after long term storage.

Description

Source	Plant-based extracts
Form	Lyophilized powder
Color	Light yellow to pale white
Solubility (Turbidity)	About 6 mg/mL in sterile water Slight cloudy in water solution

Preparation and Storage

Reconstitution	Reconstitute in sterile basal media intended to use with indicated volume (Recommend Final Conc. : 4 mg/mL)
Working Conc.	4 to 40 ug/mL of MSEP-01 as supplement in culture medium is recommended
Shipping	Shipped at ambient temperature.
Stability and Storage	2 years, below 35 °C as lyophilized powder 2 weeks, 2 to 8 °C under sterile conditions after reconstitution

Functional Profile

Promote Cell Proliferation

C2C12




Figure 1. Comparable cell growth profiles between 10% FBS and 2% FBS + MSEP-01 supplement in DMEM after 72 hours incubation in C2C12 cell (myoblast cell line).





Figure 2. Comparable cell growth profiles between 5% FBS and 2% FBS + MSEP-01 supplement in DMEM after 72 hours incubation in fibroblast culture.

Synergetic Effect with other Cell Supplement



MSEP-N: Basal medium + ITS + Albumin + Growth factor
 MSEP-C: MSEP-N + MSEP-01 (40 µg/mL)

Figure 3. MSEP has synergetic effect with ITS supplement and other growth factors to boost the cell proliferation.