Guidelines



Recombinant Collagenases (COL G and COL H) for cell isolation


1. Background

Collagen is the main structural protein of mammalian extra-cellular matrix, and it is present in different isoforms. All collagen molecules share a triple-stranded helical structure, resistant to proteases, which **can only be degraded by specific collagenases.**

The gold standard for performing cell isolation is the enzymatic digestion of the target tissue/organ by collagenases. Among collagenases the most important ones are **classes I and II** (isoforms I and II according to the nomenclature of Bond and Van Wart 1984 [1]), expressed by *C. histolyticum*, which are preferred to others, since they are **proteolytically active** practically on **all mammalian collagen isoforms**.

2. Collagenases extracted from C. histolyticum

The current technologies to produce collagenases for research and biomedical use are based on the culture of *C. histolyticum* and the subsequent protein purification thereof. The resulting 'Collagenase' is actually a **"blend"**, containing different percentage/weight ratios of the two collagenase isoforms (class I and class II), plus an ill-defined number of other proteases (clostripain, trypsin like, caseinase activity, etc.) [5].



The successful outcome of cell isolation strongly depends on activity, purity, and formulation of 'Collagenase' [4-5]. However, the extraction process of the 'Collagenase' blend from *C. hystolyticum* makes it **difficult to obtain reproducible batches**, while the presence of the proteases can contribute to a **lower stability** due to enzymatic autocatalytic processes. Low batches consistency and low stability can lead to low reproducibility and low standardization of protocols for cell isolation [6].

Collagenases

Class I

High collagenolytic activity, specifically hydrolyzing native collagen 3D-helix regions • Class II

modest activity against 3D collagen helix, acting on linear collagen regions at the motif Pro-Y-Gly-Pro [4].

Extractive Collagenases Limits in Cells Isolation

• Variability in composition of class I and II enzymes.

• Lower stability since the presence of proteases.

• **High variability in the efficiency** of cell extraction and therefore low **reproducibility** that hampers standardization of protocols.

• **Residual trypsin-like activity** despite improved purification methods.

Guidelines



3. Abiel's Recombinant Collagenases

Abiel's collagenases Class I (namely Col G) and Class II (namely Col H) of *C. histolyticum* are cloned and expressed in *E. Coli* and the tagged recombinant proteins are purified by affinity chromatography. Separate productions of the enzymes (animal-free) ensures **high yield** (>1 g/litre) and **high purity** (>99%) of COL G or COL H enzyme preparations with excellent **batch-to-batch reproducibility** [7].

In contrast to collagenase extracted from *C. histolyticum*, Abiel's recombinant products allow to prepare **customized** and **standardized mixtures** [8]. These mixtures contain a defined ratio of the two collagenases (COL G and COL H) and a defined protease. Thus, the blending is based on activity determination and not on weight ratios with ill-defined composition and activity. Taken together, the separate production allows defining the precise enzymatic activity of the final blend adapted to the use in a specific cell isolation procedure.

For cell isolation, one generic protease is required to digest the noncollagenous part of the ECM. Mixing COL G + COL H + protease, an efficient ECM digestion is obtained [9].



* either thermolysin, neutral protease/dispase, pronase. For your specific application please refer to <u>info@abielbiotech.com</u>

Benefits of Abiel's enzymes

• **Reproducible batches:** avoid batch-tobatch variation.

• **Consistency:** consistent results from each experiment. Time saving–cost saving.

• Pure collagenase (>99%): no contaminating proteases, which reduce over-digestion.

• **Reagent Flexibility:** control collagenase digestion rate (or aggressiveness) adapted to protocol requirements.

• **High stability** (up to 3 years): Lyophilised product stable at RT or cold storage; no autodegradation by proteases.

• Endotoxin-free: No off-shoot toxic effects on cells.

Applications

Human

- Langerhans' islets
- Fibroblasts (dermis regeneration)
- Mesenchymal stem cells (adipose tissue)
- Stem cells (diff. tissues)

Rat

- Langerhans' islets

Mouse

- Hepatocytes (liver regeneration)
- Mesenchymal stem cells (adipose tissue)

- Osteoblasts (bone regeneration)

Bovine

- Chondrocytes (cartilage regeneration)

References

- [1] Matsushita O. et al. (1999) *J. Bacteriol.* **181**(3): 923–933.
- [2] Philominathan S.T. et al. (2009) J. Biol. Chem. 284(16) : 10868.
- [3] Matsushita, O. Et al. (1994) J. Bacteriol.176: 149-156.
- [4] Wolters G.H. et al. (1995) *Diabetes* **44**(2) : 227-233.
- [5] Johnson P.R. et al. (1996) *Cell Transplant* **5**: 437-452.
- [6] Fermo I. et al. (2007) *Transplantation* 84(12): 1568-1575.
 [7] Salamone M. et al. (2012) *Chem. Eng. Trans.* 27: 259-264
- [8] Breite A.G. et al. (2011) *Transplant Proc.* **43**(9) : 3171-3175.
- [9] Salamone M. et al. (2014) *Chem. Eng. Trans.* **38**: 247-252.

For suggestions about your specific protocol or application of COL G and COL H, contact us:

Monica Salamone, PhD R&D Manager m.salamone @abielbiotech.com Alessandra Lo Cicero, PhD Product Manager a.locicero@abielbiotech.com

Tel. +39 091 6615633 **Fax** +39 091 661 5628