

PRODUCT INFORMATION

Anti-HMGB1 Antibody (HRP)

Catalog number	CM001-0.1EAD / 1EAD
Package	100 ul / 1 mL
Description	High-mobility group protein B1 (HMGB1) is also known as high-mobility group protein 1 (HMG-1) and amphoterin. HMGB1 is a non-histone nuclear protein with a dual function. In the nucleus, HMGB1 is a DNA binding protein that associates with chromatin and can bend DNA. Binds preferentially single-stranded DNA. Involved in V(D)J recombination by acting as a cofactor of the RAG complex. Acts by stimulating cleavage and RAG protein binding at the 23 bp spacer of conserved recombination signal sequences (RSS). Outside the cell, HMGB1 displays cytokine activity and can promote inflammation, serving as a mediator in models of shock and arthritis.
Product type	Monoclonal, Antibody pair
Host	Mouse
Reactivity	Recognizes HMGB1 protein, do not cross-react with HMGB2. Recognize HMGB1 in ELISA, when monoclonal antibody was paired with Anti-HMGB1 Antibody (cat. CM001-0.1EAC/1EAC).
Conjugation	HRP
Form	Liquid
Storage buffer	HRP stabilizer
Storage	Store at 4°C for two weeks. Store at -20°C for one year. Protect from light. Avoid freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	ELISA
Application Note	N/A
Manual	1:1000-1:5000 for ELISA. Note: Application concentration may be various determined by the end user. For Research Use Only.

For Research Use Only.